

**EvaCare Vaginal
Weights Order Form:**

Order by Phone:

866.839.9260
425.242.6373

Order by Fax:

425.497.1045

Order by Mail:

Personal Med
8531 154th Ave. NE #150
Redmond, WA 98052 USA

Questions, call:

866-839-9260

Please send the following:

Weight set(s) @ \$49.95 ea = _____ Name _____
☐ 7-10 Day Shipping = _____ Free* _____ Address _____
☐ Priority Shipping = \$10.00 _____ City/State/Zip _____
Sales Tax (WA only) = _____ Phone _____ Email _____
Order Total _____

Method of Payment (check or credit cards, no cash or HSA cards):

☐ Check payable to Personal Med
☐ Money Order or Cashier's Check
☐ Visa ☐ MasterCard ☐ American Express ☐ Discover

Card Number _____ Exp. Date _____

Signature _____ Security Code _____

*Continental US Only



Tel 866.839.9260 or +1.425.242.6373
Fax +1.425.497.1045
www.personalmed.com

U.S. Patent Nos. 5,554,092 and 6,068,581
EvaCare is a registered trademark of Personal Medical Corp
©2020-2022 Personal Medical Corp.
All rights reserved.
84411283 Rev. B January 2022



EvaCare
Vaginal Weights



*A Different Kind of
Weight Training*



EvaCare
Vaginal Weights



*Bladder control is a problem for millions of women . . .
Now you can take control!*

What Causes Bladder Control Problems?

One in four adult women experiences some type of bladder control problem, whether it's daily leakage or leakage with laughing or coughing. If you are one of these women, you are not alone.

Bladder control problems can be caused by weak muscles, specifically the pelvic floor muscles extending from the pelvic bone to the tailbone. Childbirth and aging can weaken these muscles; however, bladder control problems aren't limited to older women. Young women can also experience these problems. So how do you prevent or improve bladder control problems?



Pelvic Muscle Exercise Can Help

Many women have improved their bladder control by using pelvic muscle exercises (also called PMEs or Kegel exercises). These exercises strengthen and tone the pelvic muscles; however, they must be done correctly and consistently. EvaCare Vaginal Weights can help you do just that.

EvaCare Vaginal Weights Make Pelvic Muscle Exercise Easy

Here's how it works. EvaCare Vaginal Weights consist of four (4) graduated weights and one medical-grade plastic cone. Place one weight in the cone. Use the weight that's right for you (usually the lightest to start). Insert the weighted cone into your vagina - just like inserting a tampon. Your pelvic muscles will contract to hold it in place. To build strength and tone, keep the cone in place for up to 15 minutes at a time, twice a day. As your muscles get stronger, you can increase the weight inside the cone.

That's all there is to it. This is the easiest exercise program you'll ever do and it's virtually mistake-proof! Moving from one weight to the next allows you to track progress and literally see yourself improving. Studies have shown that pelvic muscle exercises like EvaCare Vaginal Weights can improve bladder control for many women in as little as 4-6 weeks.

Do not use vaginal weights if any of the following apply:

- You have any known or suspected vaginal, pelvic or genital diseases;
- You have a known or suspected pregnancy; or
- You have had a child or pelvic surgery in the last 6 weeks

Consult a doctor if you have an IUD, prolapse or any concerns about using vaginal weights.

Even though it's produced to medical-grade standards, the EvaCare Vaginal Weight set is very reasonably priced.



what we're putting in our bodies. Personal Med is an FDA-registered manufacturer and we produce our vaginal weights to medical-grade standards. We use biocompatible plastic, surgical (non-bacterial) suturing material and stainless steel weights, not lead like some other products.

EvaCare Vaginal Weights contain no latex.

