

eva care

# EVACARE PESSARIES





## SV GELLHORN SHORT

### *with drains*

The flexible short stem gellhorn pessary is available with drainage holes. The devices are made from medical grade silicone.

The gellhorn pessary is indicated for use for the support of a second to third-degree prolapse or procidentia.

mm	Inch	Code No.
42	1.50	SvG1.50SHD #0
45	1.75	SvG1.75SHD #1
51	2.00	SvG2.00SHD #2
57	2.25	SvG2.25SHD #3
64	2.50	SvG2.50SHD #4
70	2.75	SvG2.75SHD #5
76	3.00	SvG3.00SHD #6
83	3.25	SvG3.25SHD #7
89	3.50	SvG3.50SHD #8
95	3.75	SvG3.75SHD #9



## SV GELLHORN

### *with drains*

The flexible gellhorn pessary is available with drainage holes. The devices are made from medical grade silicone.

The gellhorn pessary is indicated for use for the support of a second to third-degree prolapse or procidentia.

mm	Inch	Code No.
42	1.50	SvG1.50D #0
45	1.75	SvG1.75D #1
51	2.00	SvG2.00D #2
57	2.25	SvG2.25D #3
64	2.50	SvG2.50D #4
70	2.75	SvG2.75D #5
76	3.00	SvG3.00D #6
83	3.25	SvG3.25D #7
89	3.50	SvG3.50D #8
95	3.75	SvG3.75D #9



## SV RING

### *with support*

The flexible SV ring pessary with support is constructed of medical grade silicone. The SV ring pessary with support is indicated for use in cases of first-degree or second-degree uterine prolapse. The SV ring pessary with support is indicated for use for uterine prolapse with accompanying cystocele.

mm	Inch	Code No.
57	2.25	SvR2.25S #2
64	2.50	SvR2.50S #3
70	2.75	SvR2.75S #4
76	3.00	SvR3.00S #5
82	3.25	SvR3.25S #6
89	3.50	SvR3.50S #7



## SV DISH

### *with support*

The flexible SV dish pessary with support is constructed of medical grade silicone. The SV dish pessary with support is indicated for stress urinary incontinence and minor degrees of prolapse.

mm	Inch	Code No.
51	2.00	SvDSH51S #0
57	2.25	SvDSH57S #1
64	2.50	SvDSH63S #3
70	2.75	SvDSH70S #4
76	3.00	SvDSH76S #5
82	3.25	SvDSH82S #6
89	3.50	SvDSH89S #8
95	3.75	SvDSH95S #9



## REGULAR RING

*with support*

The flexible ring pessary with support is constructed of medical grade silicone. The ring pessary with support is indicated for first- or second-degree uterine prolapse, accompanying cystocele.

mm	Inch	Code No.
43	1.75	R1.75S #0
51	2.00	R2.00S #1
52	2.25	R2.25S #2
62	2.50	R2.50S #3
67	2.75	R2.75S #4
74	3.00	R3.00S #5
80	3.25	R3.25S #6
86	3.50	R3.50S #7
94	3.75	R3.75S #8
98	4.00	R4.00S #9



## REGULAR RING

*without support*

The flexible ring pessary without support is made from of medical grade silicone. The ring pessary without support is indicated for use in cases of first- or second-degree uterine prolapse, and any accompanying cystocele.

mm	Inch	Code No.
43	1.75	R1.75 #0
51	2.00	R2.00 #1
55	2.25	R2.25 #2
62	2.50	R2.50 #3
64.5	2.50	R2.60 #3.5
67	2.75	R2.75 #4
70	2.75	R2.80 #4.5
74	3.00	R3.00 #5
80	3.25	R3.25 #6
86	3.50	R3.50 #7
94	3.75	R3.75 #8
98	4.00	R4.00 #9



## RING WITH KNOB

### *with support*

The ring pessary with knob, with a supporting membrane, has the same indications as the ring and may also relieve stress incontinence.

The ring with knob pessary has two notches to allow for the pessary to be easily folded.

mm	Inch	Code No.
43	1.75	RK1.75S #0
48	2.00	RK2.00S #1
56	2.25	RK2.25S #2
61	2.50	RK2.50S #3
68	2.75	RK2.75S #4
74	3.00	RK3.00S #5
81	3.25	RK3.25S #6
87	3.50	RK3.50S #7



## RING WITH KNOB

### *without support*

The ring pessary with knob, without a supporting membrane, has the same indications as the ring and may also relieve stress incontinence.

The ring with knob pessary has two notches to allow for the pessary to be easily folded.

mm	Inch	Code No.
43	1.75	RK1.75 #0
48	2.00	RK2.00 #1
56	2.25	RK2.25 #2
61	2.50	RK2.50 #3
68	2.75	RK2.75 #4
74	3.00	RK3.00 #5
81	3.25	RK3.25 #6
87	3.50	RK3.50 #7



## GELLHORN SHORT STEM

### *with drains*

The flexible short stem gellhorn pessary is available with drainage holes. The devices are made from medical grade silicone.

The gellhorn pessary is indicated for use for the support of a second to third-degree prolapse or procidentia.

mm	Inch	Code No.
42	1.50	G1.50SHD #0
45	1.75	G1.75SHD #1
51	2.00	G2.00SHD #2
57	2.25	G2.25SHD #3
64	2.50	G2.50SHD #4
70	2.75	G2.75SHD #5
76	3.00	G3.00SHD #6
83	3.25	G3.25SHD #7
89	3.50	G3.50SHD #8
95	3.75	G3.75SHD #9



## GELLHORN

### *with drains*

The flexible gellhorn pessary is available with drainage holes. The devices are made from medical grade silicone.

The gellhorn pessary is indicated for use for the support of a second to third-degree prolapse or procidentia.

mm	Inch	Code No.
42	1.50	G1.50D #0
45	1.75	G1.75D #1
51	2.00	G2.00D #2
57	2.25	G2.25D #3
64	2.50	G2.50D #4
70	2.75	G2.75D #5
76	3.00	G3.00D #6
83	3.25	G3.25D #7
89	3.50	G3.50D #8
95	3.75	G3.75D #9



## CUBE with drains

The flexible cube pessary with drains is constructed of medical grade silicone.  
The cube pessary with drains is indicated for use in the support of a third-degree prolapse, including procidentia, as well as cystocele and rectocele.

mm	Inch	Code No.
25	1	CU25D #0
29	1 1/8	CU29D #1
33	1 5/16	CU33D #2
37	1 7/16	CU37D #3
41	1 5/8	CU41D #4
45	1 3/4	CU45D #5
50	1 15/16	CU50D #6
56	2 3/16	CU56D #7
63	2 1/2	CU63D #8
70	2 3/4	CU70D #9
75	3	CU75D #10



## CUBE without drains

The flexible cube pessary without drains is constructed of medical grade silicone elastomer.  
The cube pessary without drains is indicated for use in the support of a third-degree prolapse, including procidentia, as well as cystocele and rectocele.

mm	Inch	Code No.
25	1	CU25 #0
29	1 1/8	CU29 #1
33	1 5/16	CU33 #2
37	1 7/16	CU37 #3
41	1 5/8	CU41 #4
45	1 3/4	CU45 #5
50	1 15/16	CU50 #6
56	2 3/16	CU56 #7
63	2 1/2	CU63#8
70	2 3/4	CU70 #9
75	3	CU75 #10



## DONUT

The flexible donut pessary is constructed of molded medical grade silicone elastomer. The donut pessary is indicated for use for third-degree prolapse as well as cystocele and rectocele.

mm	Inch	Code No.
50	2.00	D2.00 #0
56	2.25	D2.25 #1
64	2.50	D2.50 #2
69	2.75	D2.75 #3
75	3.00	D3.00 #4
82	3.25	D3.25 #5
87	3.50	D3.50 #6
93	3.75	D3.75 #7



## SHAATZ

*with drains*

The flexible shaatz pessary is constructed of molded medical grade silicone elastomer. The shaatz pessary is indicated for use for the support of first or second-degree prolapse and an accompanying cystocele.

mm	Inch	Code No.
38	1.50	SH1.50 #0
45	1.75	SH1.75 #1
51	2.00	SH2.00 #2
57	2.25	SH2.25 #3
64	2.50	SH2.50 #4
70	2.75	SH2.75 #5
76	3.00	SH3.00 #6
81	3.25	SH3.25 #7
87	3.50	SH3.50 #8



## DISH

*with support*

The flexible dish pessary with support is made of medical grade silicone. The dish pessary with support is indicated for stress urinary incontinence and minor degrees of prolapse.

mm	Inch	Code No.
50	1.75	DSH50S #0
55	2.00	DSH55S #1
60	2.25	DSH60S #2
65	2.50	DSH65S #3
70	2.75	DSH70S #4
75	3.00	DSH75S #5
80	3.25	DSH80S #6
85	3.50	DSH85S #7



## DISH

*without support*

The flexible dish pessary without support is made of medical grade silicone. The dish pessary without support is indicated for stress urinary incontinence and minor degrees of prolapse.

mm	Inch	Code No.
50	1.75	DSH50 #0
55	2.00	DSH55 #1
60	2.25	DSH60 #2
65	2.50	DSH65 #3
70	2.75	DSH70 #4
75	3.00	DSH75 #5
80	3.25	DSH80 #6
85	3.50	DSH85 #7



## CUP

*with support*

The cup pessary, with a supportive membrane, is used for a first or a mild second-degree uterine prolapse.

mm	Inch	Code No.
54	2.00	CP55S #1
58	2.25	CP60S #2
63	2.50	CP65S #3
68	2.75	CP70S #4
75	3.00	CP75S #5
79	3.25	CP80S #6



## CUP

*without support*

The cup pessary, without a supportive membrane, is used for a first or a mild second-degree uterine prolapse.

mm	Inch	Code No.
54	2.00	CP55 #1
58	2.25	CP60 #2
63	2.50	CP65 #3
68	2.75	CP70 #4
75	3.00	CP75 #5
79	3.25	CP80 #6



## OVAL

*with support*

The flexible oval pessary is made from medical grade silicone. The oval pessary is indicated for the support of first or second-degree prolapse as well as an accompanying cystocele. It is specifically designed to fit a narrow vaginal vault.

mm	Inch	Code No.
51	2.00	OV2.00S #1
57	2.25	OV2.25S #2
64	2.50	OV2.50 #3
70	2.75	OV2.75 #4
75	3.00	OV3.00 #5
82	3.25	OV3.25 #6
89	3.50	OV3.50 #7
95	3.75	OV3.75 #8
101	4.00	OV4.00 #9



## GEHRUNG

*with support*

The gehrung pessary with support is made from medical grade silicone surrounding a flexible stainless steel wire, which allows manual shaping. The gehrung pessary with support is indicated for use for the support of cystoceles and rectoceles, as well as second to third-degree uterine prolapse.

mm	Inch	Code No.
50	1.75	GH50S #0
55	2.00	GH55S #1
60	2.25	GH60S #2
65	2.50	GH65S #3
70	2.75	GH70S #4
75	3.00	GH75S #5
80	3.25	GH80S #6
85	3.50	GH85S #7
90	3.75	GH90S #8
95	4.00	GH95S #9



## MAR-LAND

### *with support*

The flexible mar-land pessary with support is made from medical grade silicone elastomer. The mar-land pessary is indicated for stress urinary incontinence and minor prolapse.

mm	Inch	Code No.
57	2.25	M2.25S #2
64	2.50	M2.50S #3
70	2.75	M2.75S #4
76	3.00	M3.00S #5
83	3.25	M3.25S #6
89	3.50	M3.50S #7
95	3.75	M3.75S #8



## MAR-LAND

### *without support*

The flexible mar-land pessary without support is made from medical grade silicone. The mar-land pessary is indicated for stress urinary incontinence and minor prolapse.

mm	Inch	Code No.
57	2.25	M2.25 #2
64	2.50	M2.50 #3
70	2.75	M2.75 #4
72	3.00	M3.00 #5
83	3.25	M3.25 #6
89	3.50	M3.50 #7
95	3.75	M3.75 #8



# HODGE

*with support*

The folding hodge pessary with support is made from medical grade silicone surrounding a flexible stainless steel wire, to allow manual shaping. The hodge pessary is indicated for a first to second-degree prolapse, a cystocele, stress urinary incontinence and an incompetent cervix or uterine retroversion.

mm	Inch	Code No.
65	2.50	HD65S #0
70	2.75	HD70S #1
75	3.00	HD75S #2
80	3.25	HD80S #3
85	3.50	HD85S #4
90	3.75	HD90S #5
95	4.00	HD95S #6
100	4.25	HD100S #7
105	4.50	HD105S #8
110	4.75	HD110S #9



# HODGE

*without support*

The folding hodge pessary without support is made from medical grade silicone surrounding a flexible stainless steel wire, to allow manual shaping. The hodge pessary is indicated for a first to second-degree prolapse, a cystocele, stress urinary incontinence and an incompetent cervix or uterine retroversion.

mm	Inch	Code No.
65	2.50	HD65 #0
70	2.75	HD70 #1
75	3.00	HD75 #2
80	3.25	HD80 #3
85	3.50	HD85 #4
90	3.75	HD90 #5
95	4.00	HD95 #6
100	4.25	HD100 #7
105	4.50	HD105 #8
110	4.75	HD110 #9

eva care



## NOT ALL PESSARIES ARE ALIKE ...

### *Feel the difference*

Personal Med sells the 10 or so different types of pessaries, but probably the two most important are the Rings and Dishes:

- Rings are used for milder prolapse (first to second degree) and are by far the most popular kind of pessary.
- Dishes are used when milder prolapse also results in urinary incontinence, which is unfortunately quite common.

Because of their popularity and clinical utility, Personal Med offers a special version (SV) of these two types of pessaries. Our SV pessaries, which are recognizable from their modified Yin symbol, have a different silicone formulation that is softer than other pessaries. Many clinicians and patients report that EvaCare SV pessaries are easier to insert and remove.

Try them and see how you feel.


Personal Med's mission is to advance the pelvic health of women with incontinence and /or prolapse. We do this by providing best-in-class medical devices supported by a knowledgeable, caring staff.

# EVACARE SIZING CHARTS

Please call or email to receive any of these useful aids for your clinical practice.

Phone: (425) 242 6373

Email: [customerservice@personalmed.com](mailto:customerservice@personalmed.com)



The logo for 'eva care' features the text 'eva care' in a purple, lowercase, sans-serif font. The text is centered within a white, teardrop-shaped area that is partially enclosed by a teal-colored shape above and a purple-colored shape below, creating a stylized, organic frame.

eva care

Personal Med Corp  
8531 154th Ave NE  
REDMOND WA 98052  
(425) 242 6373

The logo for 'personalmed' features a stylized icon of three white, teardrop-shaped elements arranged in a slight arc above the text.

**personalmed**  
*advancing pelvic health*<sup>®</sup>